

Pipettes



A3005



A3008



A3021

A3018

Pipette Pump

| Cat No. | Description | Packing Unit |
|---------|-------------|---------------------|
| A3005 | 2ml | Independent packing |
| A3006 | 10ml | Independent packing |
| A3007 | 25ml | Independent packing |

Lab Adjustable-Volume Pipettes

- Fully autoclavable—entire pipette can be steam autoclaved.
- Extremely low plunger forces.
- Soft-touch tip ejection.
- Ergonomic comfortable lightweight handle.
- Super blow-out piston.

| Cat No. | Volume Range(µl) | Increment(µl) | Test Volume (µl) | Inacc.(%) | Impc.(%) | Tip |
|---------|------------------|---------------|------------------|-----------|----------|----------------|
| A3008 | 0.5-10µl | 0.1 | 10µl | ±3.0 | ≤2.0 | 10µl |
| | | | 5µl | ±6.0 | ≤3.0 | |
| | | | 1µl | ±10.0 | ≤3.0 | |
| A3009 | 5-50µl | 0.5 | 50µl | ±2.4 | ≤1.5 | 300µl 350µl |
| | | | 25µl | ±2.8 | ≤2.0 | |
| | | | 5µl | ±6.0 | ≤3.0 | |
| A3010 | 10-100µl | 1 | 100µl | ±1.5 | ≤1.0 | 300µl 350µl |
| | | | 50µl | ±2.2 | ≤1.5 | |
| | | | 10µl | ±3.0 | ≤2.0 | |
| A3011 | 100-1000µl | 5 | 1000µl | ±1.2 | ≤0.5 | 1000µl |
| | | | 500µl | ±1.5 | ≤0.5 | |
| | | | 100µl | ±1.5 | ≤1.0 | |
| A3012 | 1000-5000µl | 50 | 5000µl | ±1.0 | ≤0.5 | 5000µl |
| | | | 2500µl | ±1.0 | ≤0.5 | |
| | | | 1000µl | ±1.2 | ≤0.5 | |

Lab Fixed Volume Pipette

| Cat No | Volume(µl) | Inacc. % | Impc. % |
|--------|------------|----------|---------|
| A3013 | 5 | ±4.8 | ≤3.0 |
| A3014 | 10 | ±3.0 | ≤2.0 |
| A3015 | 100 | ±1.5 | ≤0.8 |
| A3016 | 500 | ±1.2 | ≤0.7 |
| A3017 | 1000 | ±1.0 | ≤0.5 |

Lab Multichannel Pipettes

- 8-and 12-channel pipettes with volume range of 0.5 up to 300µl.
- Fully autoclavable—entire pipette can be steam autoclaved.
- Extremely low plunger forces.
- Soft-touch tip ejection.
- Ergonomic comfortable lightweight handle.
- Optional fit with Sureone pipette tips.
- Super blow out piston.

| Cat No | Volume Range | Increments | Each |
|-------------------|--------------|------------|---------------------|
| 8-Channel | | | |
| A3018 | 0.5-10µl | 0.1µl | Independent packing |
| A3019 | 5-50µl | 0.5µl | Independent packing |
| A3020 | 50-300µl | 5µl | Independent packing |
| 12-Channel | | | |
| A3021 | 0.5-10µl | 0.1µl | Independent packing |
| A3022 | 5-50µl | 0.5µl | Independent packing |
| A3023 | 50-300µl | 5µl | Independent packing |