

Autoclave

Performance :

- > Microprocessor controller and temperature calibration function allows for high precision temperature adjustment.
- > Sterilization temperature at 134°C
- > Automatic sterilization system for unattended operation.



BA-18

Model		BA-8	BA-12	BA-16	BA-18	BA-23
Cat No.	110V 60HZ	70651000	70651200	70651400	70651600	70651800
	220V 50HZ	70652000	70652200	70652400	70652600	70652800
Capacity		8 L	12 L	16 L	18 L	23 L
Consumed Power		950 W	1600 W	1600 W	1700 W	1900 W
Working Pressure		0.22Mpa				
Working Temperature		134°C				
Chamber Volume		φ170×320(mm)	φ 200×360(mm)	φ 200×360(mm)	φ 247×350(mm)	φ 247×450(mm)
Weight		32kg	43kg	44kg	48kg	53kg
		71LB	95LB	97LB	107LB	118LB
Dimension(mm/inch)		640×480×430	690×560×490	690×560×490	690×590×530	790×590×540
		25×19×17"	27×22×19"	27×22×19"	27×23×21"	31×23×21"

Convenience :

- Fully automatic computer control, “LCD” dynamic display, fault self - diagnosis, display the error code;
- Preset sterilization procedures, and adjustable drying time;
- SUS304 stainless steel inner and the sealing cover ,double pressure protection safety door lock;
- Built-in independent steam generator, fast saturated steam, to insure the high temperature steam, reaches every corner;
- Simple and convenient operating interface, LCD displays the status information in real time;
- Open type top water tank, filling and cleaning is convenient for user.

Safety :

- Machine runs ONLY when doors are closed & locked;
- Safety pressure protection lock, Automatically, inducting the pressure inside the sterilizing cavity, can't open door during sterilization process;
- With liquid sensor for water level, to remind user to add water;



BA-18V

Model		BA-8V	BA-12V	BA-16V	BA-18V	BA-23V
Cat No.	110V 60HZ	70651100	70651300	70651500	70651700	70651900
	220V 50HZ	70652100	70652300	70652500	70652700	70652900
Capacity		8 L	12 L	16 L	18 L	23L
Consumed Power		950 W	1750 W	1750 W	1800 W	2000 W
Working Pressure		0.22Mpa				
Working Temperature		134℃				
Chamber Volume		φ170×320(mm)	φ 200×360(mm)	φ 230×360(mm)	φ 247×350(mm)	φ 247×450(mm)
Weight		37kg	46kg	47kg	52kg	58kg
		82LB	102LB	104LB	116LB	129LB
Dimension(mm/inch)		640×480×430	690×560×490	690×560×490	690×590×530	790×590×540
		25×19×17″	27×22×19″	27×22*19″	27×23×21″	31×23×21″

Autoclave

Digital Display – Vertical Automatic Autoclave

- Fully automatic control system with memory function.
- 2 programs preset, Automatic control function stronger.
- Self-inflating type sealing, safe and reliable.
- Safety interlocking pressure system, To avoid wrong operation cause by steam burns.
- Four protection functions: self discharge over pressure, over pressure cut off the heating, anti dry burning, low water level alarm.
- Microprocessor fully automatic control the discharge, emissions and frequency can be regulated freely; sterilization degree and cycle can be self balance.
- With manual exhaust function, solve the problem of negative pressure cover opening difficult and special needs; (such as slow exhaust).
- Automatically shut off with beep reminding after sterilization.
- With standard test interface.
- The main parts are acid and alkali resistant, corrosion resistant and high quality SUS304 stainless steel materials.
- Patent technology: easy clean up the water level control device, more adaptable to hard water (tap water).
- Optional automatic water (connected with the pressure water, such as tap water).
- With optional external printer.



VA-75B

Model		VA-30B	VA-40B	VA-50B	VA-75B	VA-100B
Cat No.	110V 60HZ	70661100	70661200	70661300	70661400	70661500
	220V 50HZ	70662100	70662200	70662300	70662400	70662500
Capacity		30 L	40 L	50 L	75L	100L
Consumed Power		2.0 KW	2.0 KW	3.5 KW	3.5KW	3.5KW
Working Pressure		0.22Mpa				
Working Temperature		134 °C				
Chamber Volume		φ300×400	φ300×600	φ400×525	φ400×625	φ400×725
Basket Dimension		φ280×360×1	φ280×280×2	φ360×240×2	φ360×280×2	φ360×240×3
Weight		53kg	55kg	58kg	77kg	120kg
		117 LB	121LB	128LB	170LB	265LB
Dimension(mm/inch)		410×410×850	410×410×1050	500×500×1080	500×500×1150	500×500×1280
		16×16×33"	16×16×41"	20×20×43"	20×20×45"	20×20×50"

Portable steam sterilizer series

PS-18A/A+model

- All made of SUS304 stainless steel material, durable.
- Double scale indication pressure gauges, clear to read.
- Self-inflating type sealing, safe and reliable.
- Auto release at overpressure.
- Simple operation, safe and reliable.
- Optional manual knob timing 60min.



PS-18B/B+ Automation type

- All made of SUS304 stainless steel material, durable.
- Microcomputer control, digital display, sterilization temperature and time can be adjusted.
- Auto release at overpressure.
- Double scale indication pressure gauges, clear to read.
- Self-inflating type seal, safe and reliable.
- Simple operation, safe and reliable.



PS-18 C/C+ (Portable) Electric or LPG heated

- All made of SUS304 stainless steel material, durable.
- Automatic control the pressure and temperature of sterilization.
- Auto release at overpressure.
- Double scale indication pressure gauges, clear to read.
- Self-inflating type seal, safe and reliable.
- Electric heating and external heating both can use.
- Simple operation, safe and reliable.



Model		PS-18A	PS-18A+	PS-18B	PS-18B+	PS-18C	PS-18C+
Cat No.	110V 60HZ	70671100	70671200	70671300	70671400	70671500	70671600
	220V 50HZ	70672100	70672200	70672300	70672400	70672500	70672600
Volume		18L	24L	18L	24L	18L	24L
Consumed Power		2000W	2000W	2000W	2000W	2000W	2000W
Working Pressure		0.165 Mpa					
Working Temperature		126℃					
Sterilizing room dimensions		Φ280×235	Φ280×340	Φ280×235	Φ280×340	Φ280×235	Φ280×340
Packing size(mm/inch)		410×410×550	425×425×660	410×410×550	425×425×660	400×400×440	400×400×540
		16×16×22"	17×17×26"	16×16×22"	17×17×26"	16×16×17"	16×16×21"
Weight		14.5kg	19kg	15.5kg	22.5kg	12kg	17kg
		32LB	42LB	34LB	50LB	27LB	37LB