

Biological Microscope



AMB-17



AMB-17C

PARAMETERS	BIOLOGICAL MICROSCOPE	AMB-17	AMB-17C
Total Magnification	40X-1600X	●	●
Eyepieces	WF10X WF16X	●	●
Objectives	4X,10X,40XS,100XS	●	●
Eyepieces Tube	Binocular 45° inclined Tube Rotatable 360°	●	●
Illumination	6V/20W Halogen Lamp	●	●
Coarse Adjustment Range	30mm	●	●
Fine Adjustment Range	30mm	●	●
Working Stage	Double Layer Mechanical Stage:130*140mm Moving Range:75mm*45mm	●	●
Condenser	NA1,25 Abee Condenser with Changeable Diaphragm,Moving Scale And Filter	●	●
Package Size (cm/inch)	32×35×49 12.6×13.8×19.3"	●	●
Weight	9kg / 20 LB	●	●
Camera			●
Cat No.	100V-240V	70800100	70800200

Biological Microscope



AMB-01



AMB-02

Parameters	Biological Microscope	AMB-02 (Binocular)	AMB-01 (Monocular)
Total Magnification	40X-1600X	●	●
Eyepieces	WF10X WF16X	●	●
Objectives	4X,10X,40X(S),100X,(S,Oil)	●	●
Eyepieces Tube	Binocular 30° inclined Tube Rotatable 360° Monocular 45° inclined Tube Rotatable 360°	●	●
Illumination	6V/20W Halogen Lamp	●	●
Coarse Adjustment Range	20mm/14mm	20mm	14mm
Working Stage	Double Layer Mechanical Stage:125×130mm Moving Range:60mm×30mm	●	●
	Double Layer Mechanical Stage:110×125mm Moving Range:60mm×20mm		
Condenser	NA1,25 Abee Condenser with Changeable Diaphragm,Moving Scale And Filter	●	●
Package Size (cm/inch)	50×34×50 /19.7×13.4×19.7" 36×26×21 /14.2×10.2×8.3"	●	●
Weight	12kg/26.4LB 3.5kg/7.7 LB	●	●
Cat No./ Power	100-240V	70800300	70800400